



# DiViDu-Next

## **Installation guide**


### **Windows**

## Instructions to download and run the DiViDu-Next files on Windows:


1. The latest versions of the code files are available on GoogleDrive at the following link:

[https://drive.google.com/file/d/1A8kg3Ktv9qi6hLaDsx5A0ALAo52yQEkf/view?usp=drive\\_link](https://drive.google.com/file/d/1A8kg3Ktv9qi6hLaDsx5A0ALAo52yQEkf/view?usp=drive_link)


2. This link will take you to a download page.
3. Here click on the blue download button.



4. Wait for the dividunext.zip file to download.
5. Once the folder is downloaded, locate the "DividuNext.zip" file on your computer.
6. Right click on the "DividuNext.zip" file and select "Extract All".




7. Select “Extract” in the window that pops up.




8. Wait for the files to be extracted. This may take a minute or two.


9. Navigate to and open the “DividuNext” folder.
10. Inside the main folder, in the search bar located in top right corner type “main.exe”.




11. Now you will see the “main.exe” application file. Double click on this application to launch the experiment window.
12. You may this screen before the experiment window opens. No need to worry. Windows simply does not recognize this new application. Click on “More info”.



13. Next click on “Run anyway” at the bottom.



14. Now you should see the experiment window and you are ready to start.



15. If you close the Dividu-next application window you can open the application again by double clicking on the “main.exe” application file.
16. The experiment outputs will be stored in the “reports” folder inside the “DividuNext” folder. You can type “reports” in the top right search bar to easily locate the reports.